

System Name	Product Name	Use with main product as combination (Product or Accessory)					Instrument Set code	Instrument Set Name	Instrument Set Details	Mention if used on joints (Implants used on joints will affect the adjacent body part)
	Plate Length	Screws With Quantity	Pins (If Applicable) With Quantity	Wire With Quantity	Bolt (If Applicable) With Quantity	Fixators (If Applicable) With Quantity				
3.5mm Antrolateral Distal Tibia Locking – Left 5 hole 7 hole 9 hole 11 hole 13 hole 15 hole	82mm 108mm 134mm 160mm 186mm 212mm	3.5mm Locking Screw (10mm – 40mm) 3.5mm Cortical Screw (10mm – 40mm) 4.0mm Locking Screw (20mm-70mm)	-	K-wire Double Trocar 2.0mm	-	-	BO.511.SET	Antrolateral Distal Tibia Locking Instrument Set	Lcp Drill Sleeve 3.2mm & 2.7mm, Drill Bit 3.2mm & 2.7mm, 2.5mm Tip Screw Driver, Bone Tap 3.5mm, Depth Gauge, Torque Screw Driver 3.5mm	No Effect on Bone & Adjacent Body Part