

| System Name   | Product Name                     | Use with main product as combination (Product or Accessory)             |                                    |                            |                                    |  | Instrument Set code | Instrument Set Name              | Instrument Set Details  | Mention if used on joints (Implants used on joints will affect the adjacent body part) |
|---|----------------------------------|---|------------------------------------|----------------------------|------------------------------------|--|---------------------|----------------------------------|---|--|
|   | Plate                            | Screws With Quantity  | Pins (If Applicable) With Quantity | Wire With Quantity         | Bolt (If Applicable) With Quantity | Fixators (If Applicable) With Quantity |                     |                                  |   |  |
| 3.5mm Olecranon Locking Left<br>4 hole<br>6 hole<br>8 hole<br>10 hole | 126mm<br>152mm<br>178mm<br>204mm | 3.5mm Locking Screw (10mm – 40mm)<br>3.5mm Cortical Screw (10mm – 40mm) | -                                  | 2.0mm K-wire Double Trocar | -                                  | -                                      | BO.509.SET          | Olecranon Locking Instrument Set | Lcp Drill Sleeve<br>3.2mm & 2.7mm, Drill Bit<br>3.2mm & 2.7mm,<br>2.5mm Tip Screw Driver,<br>Bone Tap<br>3.5mm,<br>Depth Gauge,<br>Torque Screw Driver<br>3.5mm | No Effect on Bone & Adjacent Body Part   |