

System Name	Product Name	Use with main product as combination (Product or Accessory)					Instrument Set code	Instrument Set Name	Instrument Set Details	Mention if used on joints (Implants used on joints will affect the adjacent body part)
	Plate	Screws With Quantity	Pins (If Applicable) With Quantity	Wire With Quantity	Bolt (If Applicable) With Quantity	Fixators (If Applicable) With Quantity				
3.5mm Posteromedial Proximal Tibia Locking – Right 3 hole 4 hole 5 hole 6 hole 7 hole 8 hole 9 hole 10 hole 11 hole 12 hole	83mm 96mm 109mm 122mm 135mm 148mm 161mm 174mm 187mm 200mm	3.5mm Locking Screw (10mm – 40mm) 4.0mm Locking Screw (20mm – 85mm) 3.5mm Cortical Screw (10mm – 40mm)	-	1.5mm K-Wire Double Trocar	-	-	BO.273.SET	Posteromedial Proximal Tibia Locking Instrument Set	Lcp Drill Sleeve 3.2mm & 2.7mm, Drill Bit 3.2mm & 2.7mm, 2.5mm Tip Screw Driver, Bone Tap 3.5mm, Depth Gauge, Torque Screw Driver 3.5mm	No Effect on Bone & Adjacent Body Part