

| System Name                             | Product Name                                 | Use with main product as combination (Product or Accessory) |                                    |                    |                                     |  | Instrument Set code | Instrument Set Name    | Instrument Set Details  | Mention if used on joints (Implants used on joints will affect the adjacent body part) |
|---|--|---|------------------------------------|--------------------|-------------------------------------|--|---------------------|------------------------|---|--|
|   | Nail   | Screws With Quantity  | Pins (If Applicable) With Quantity | Wire With Quantity | Bolt (If Applicable) With Quantity  | Fixators (If Applicable) With Quantity |                     |                        |   |  |
| Long Proximal Femur Nail Right 135° 9mm | 34cm<br>36cm<br>38cm<br>40cm<br>42cm<br>44cm | -   | -                                  | Guide Wire 2.5mm   | 4.9mm Interlocking Bolt (20mm-60mm) | -                                      | BO.712              | PFN-TFN Instrument Set | PFN-TFN Zig, Sleeve, Drill Bit, Conical Bolt, Reamer, Bone Awl, Tissue Protector, Extractor Rod, Handle, Hammer, Depth Gauge, Spanner, Screw Driver, Cap. | NO   |
| 10mm                                    | 34cm<br>36cm<br>38cm<br>40cm<br>42cm<br>44cm |   |                                    |                    | 8.0mm PFN Bolt (40mm-120mm)         |  |                     |                        |   |  |
| 11mm                                    | 34cm<br>36cm<br>38cm<br>40cm<br>42cm<br>44cm |   |                                    |                    | 6.4mm PFN Bolt (35mm-120mm)         |  |                     |                        |   |  |

|      |  |  |  |  |  |  |  |  |  |  |
|------|--|--|--|--|--|--|--|--|--|--|
| 12mm | 34cm<br>36cm<br>38cm<br>40cm<br>42cm<br>44cm |  |  |  |  |  |  |  |  |  |
|------|--|--|--|--|--|--|--|--|--|--|

