

System Name	Product Name	Use with main product as combination (Product or Accessory)					Instrument Set code	Instrument Set Name	Instrument Set Details	Mention if used on joints (Implants used on joints will affect the adjacent body part)
	Plate	Screws With Quantity	Pins (If Applicable) With Quantity	Wire With Quantity	Bolt (If Applicable) With Quantity	Fixators (If Applicable) With Quantity				
Proximal Humerus Locking (Philos) 2 hole 3 hole 4 hole 5 hole 6 hole 7 hole 8 hole 9 hole 10 hole 11 hole 12 hole	90mm 108mm 126mm 144mm 162mm 180mm 198mm 216mm 234mm 252mm 270mm	3.5mm Locking Screw (10mm – 40mm) 4.0mm Locking Screw (20mm – 70mm) 3.5mm Cortical Screw (10mm – 40mm)	-	1.5mm K-Wire Double Trocar	-	-	BOI.156.SET	Proximal Humerus Locking (Philos) Instrument Set	Lcp Drill Sleeve 3.2mm & 2.7mm, Drill Bit 3.2mm & 2.7mm, 2.5mm Tip Screw Driver, Bone Tap 3.5mm, Depth Gauge, Torque Screw Driver 3.5mm	No Effect on Bone & Adjacent Body Part